## **ACCU-CHEK** <sup>®</sup> Guide



SIMPLY MORE FOR YOUR PATIENTS

The Accu-Chek Guide meter does more to help patients manage diabetes. It offers a range of convenient, practical features for a better testing experience.

In a survey of people with diabetes who evaluated the Accu-Chek Guide SmartPack vial and test strip:

of patients found it was easier to remove a strip compared to a standard round vial.

of patients found it was easier and faster to correctly place a small amount of blood on the end of the strip.<sup>1</sup>





## **ACCU-CHEK®** Guide



## **EASY**TESTING

**Wide dosing area:** Apply a blood drop anywhere along the entire width of the yellow edge.

Spill resistant vial: Easy to remove strip from the vial.1

**Illuminated test strip slot:** lets patients easily test day or night.

Strip ejector: for quick, easy hygienic disposal.

More insights: Easy to review your averages and results such as how often values are in or out of range.

## Accu-Chek FastClix lancing device:

1-Click lancing device with a drum for virtually pain free lancing.<sup>2</sup>



CLEAR AND SIMPLE TRACKING

The Accu-Chek Guide system can be connected to the mySugr app via *Bluetooth®* to log blood glucose values automatically and create a seamless and easy diabetes management experience.



<sup>2</sup>Kocher S, et al. Comparison of lancing devices for self-monitoring of blood glucose regarding lancing pain. J Diabetes Sci Technol. 2009;3(5):1136-1143.

